




Multipurpose Ports  
 Overcurrent Protection  
 Optional Accessories  
 AC Power Supply

2in Access Port  
 (1) Each Side

Humidifier Water Supply\*  
 Compressed Dry Air Supply\*



304 Stainless Steel Wire Shelf  
 (6) Included

Notes:

- Crated Weight: 1350lbs.
- Shelves adjustable in 1/2" increments. Minimum 2" spacing required.
- Weight Capacity: 150lbs per shelf.
- Dimensions and components are subject to change without notice.
- Some options may affect chamber performance.