

## BUILT-IN UNDERCOUNTER FREEZER

### Series Overview

The built-in undercounter freezer features a microprocessor temperature controller and display, this allows for precise temperature control, verification and recovery versus mechanical style thermostats. Models are designed with "front-breathing" ventilation which eliminates the need for clearance around the top and sides of the unit. \*Periodic manual defrosting is required

### BIUF-0420SS

One stainless steel swing door chamber designed for maximum capacity while still allowing passage through most commercial doorways.

### Performance

Unit comes standard with a temperature range of -15° C to -25° C Manual defrost. Products designed to operate in ambient room temperatures not to exceed 75° F / 50% RH.

### Control System

All series models come standard with a digital microprocessor temperature controller providing a readout of actual and set-point values. This chamber comes with audible/visual alarms, high and low temperature alarms with remote alarm contacts.

### Warranties

These chambers come with a two year parts and labor warranty.



The Premier Built-in Undercounter Refrigerators and Freezers take a large step forward in performance from past general purpose units; offering a digital temperature display and alarm module coupled with a microprocessor temperature controller.



### Construction

Freezer is constructed of stainless steel exterior.

- 2 fixed shelves
- one stainless steel swing door, right hinged
- Magnetic door gasket for positive seal with keyed door lock and recessed door handle
- Leveling legs
- Probe access port (3/8")
- High density urethane foam cabinet and door insulation
- HFC-free refrigerant (R600a), foam insulation and packaging
- Exterior dimensions: 23.75" W x 24.5"D x 33.38"H
- Shipping weight: 131 lbs
- Voltage: 115 / 60 / 1.5 Amps, 1/8 HP
- Energy consumption: 0.97 kWh/day
- UL/C-UL listed
- Pyxis®, Omnicell® and AcuDose RX® compatible
- Energy Star Certified



2945 Washington Ave. St. Louis, MO 63103

314-534-3111 | [sales@darwinchambers.com](mailto:sales@darwinchambers.com)

[www.darwinchambers.com](http://www.darwinchambers.com)

---

## Built-In Undercounter Freezer Technical Drawings