

DARWIN CHAMBERS: CHILLIE

Darwin Chambers Portable Refrigerated Incubator - the "Chillie" - is supplied with all the accessories you need for refrigeration storage almost anywhere. The Chillie is an exceptional portable cooler / freezer for vehicle use whether it is being used for out of the lab experiments, research, or lab work. The Chillie is great for bigger samples and research projects.

Construction

- Top-opening, hinged and latch lid, with durable side handles for easy carrying and transportation.
- Durable powder coated steel exterior construction.
- Interior is easy to clean, seamless plastic
- Insulation is foamed in place polyurethane, with a high r-value/inch.
- Unit completely pre-wired at factory, and includes a 12VDC car adapter and AC/DC power supply cord.
- Microprocessor LED temperature control with display
- Visual Alarms, absolute hi/lo, 0.1°C increments
- Compressor Troubleshoot Fault LED control display

Capacity

- Temperature Range: -20°C to 10°C
- Inside Dimensions: 10.5" L x 14.5" D x 10.5" W of usable space
- Weight: 53 lbs.
- Voltage Power AC (110V - 65W / 2.5A) or DC (12VC - 4.5A Car Lighter Socket)



Features

- Latch Lid
- Insulated Lid
- Tough and solid outer casing
- Functions even when tilted 30 degrees
- 1 Wire Basket
- Side Handle
- DC Power Cord
- AC Power Cord
- Air Vents
- Drain Plug
- PT100 RTD Sensor
- 2 Circulating Fans

Warranty

Darwin Chambers Company offers a 6 month warranty on all parts and labor for the Chillie.