

## DARWIN CHAMBERS: CHILLIE

Darwin Chambers Portable Refrigerated Incubator - the "Chillie" - is supplied with all the accessories you need for refrigeration storage almost anywhere. The Chillie is an exceptional portable cooler / freezer for vehicle use whether it is being used for out of the lab experiments, research, or in the lab work. The Chillie is great for bigger samples and research projects.

### Construction

- Top-opening, hinged and latch lid, with durable side handles for easy carrying and transportation.
- Durable powder coated steel exterior construction.
- Interior is easy to clean, seamless plastic
- Insulation is foamed in place polyurethane, with a high r-value/inch.
- Unit completely pre-wired at factory, and includes a 12VDC car adapter and AC/DC power supply cord.
- Microprocessor LED temperature control with display
- Visual Alarms, absolute hi/lo, 0.1°C increments
- Compressor Troubleshoot Fault LED control display

### Capacity

- Temperature Range: -20°C to 10°C
- Inside Dimensions: 11"W x 11.25"L x 14.75"H of usable space
- Weight: 53 lbs.
- Voltage Power AC (110V - 65W / 2.5A) or DC (12VC - 4.5A Car Lighter Socket)



*Catalog Number: Chillie*



### Features

- Latch Lid
- Insulated Lid
- Tough and solid outer casing
- Functions even when tilted 30 degrees
- 1 Wire Basket
- Side Handle
- DC Power Cord
- AC Power Cord
- Air Vents
- Drain Plug
- PT100 RTD Sensor
- 2 Circulating Fans

### Warranty

Darwin Chambers Company offers a 6 month warranty on all parts and labor for the Chillie.